




AQUA
SCIENCE

MILIEUMONITORING MET BIOASSAYS

Praktisch aan de slag met Microbiotests Toxkits & BioLight

- 
- The background image is an aerial photograph of a wastewater treatment plant. It features several large circular tanks, some with blue water and others with brown sludge. There are also rectangular aeration basins with metal walkways and railings. A small building with a white roof is visible in the upper left. The overall scene is industrial and green.
- The number of chemicals encountered in water, soil and sediment is increasing worldwide. (> 350.000 compounds)
 - Chemical measurements provide information about the presence of chemicals in water, when it is not yet clear whether these might be harmful to human and environmental health.
 - Effect-based monitoring with bioassays answers this concern.

MARINE DUMPING



Industrial Waste



AGRICULTURAL RUNOFF



RADIOACTIVE CONTAMINATION



EXAMPLES OF WATER POLLUTION

HOUSEHOLD WASTE



THERMAL POLLUTION



OIL SPILLS



MINING ACTIVITIES



PLASTIC WASTE




SEWAGE AND WASTEWATER



CHEMICAL
ANALYSIS



BIOASSAYS



KNOWN

UNKNOWN

BIOASSAYS

MIXTURES



BIOASSAYS

- Bioassays make use of organisms, cells and bacteria, and are an effective instrument in investigating the effects of chemicals on humans and the environment (ecotoxicity).
- Bioassays measure the combined effect of chemicals present in chemical mixtures, without the need of having information about which chemicals are present and in which concentrations.
- Bioassays can be conducted, among others, on drinking water, wastewater, surface water and groundwater.
- In vivo and in vitro bioassays




Principle of a laboratory bioassay



Observed effects : mortality, growth, reproduction, DNA damage, ...



After a given time, we count the number of individuals affected in each container



IN VIVO BIOASSAYS =
need of living organisms



Burden of Culturing Test Organisms

- **Infrastructure-intensive :**

Requires lab space, aquaria, lighting and constant temperature

- **Resource-demanding :**

Daily feeding, water changes and breeding maintenance


- **Expertise-dependent :**

Requires trained personnel and experience

- **Risk of variability :**

Sensitivity and health of organisms can vary with culturing conditions









Microbiotests TOXKITS & TEST BIOTA

**SOLUTIONS FOR
ENVIRONMENTAL TOXICITY TESTING
AND
RESEARCH**

Introduction to Microbiotests and Toxkits

- Microbiotests Toxkits : Ready-to-use bioassays starting from dormant life stages
- Organisms are hatched on-demand from cysts, resting eggs or immobilized stages
- Long-term storage, culture-independent, and easy to use
- Provide reproducible and comparable results across labs :
standardized test procedures (ISO and OECD)
- Genetic uniform and qualitative biological materials



Application areas

- Registration of new chemical substances
(detergents, cosmetics, phytosanitary products, pesticides, antifoulings, ...)
- Monitoring of wastewater from urban wastewater treatment plants and industrial effluents (WWTP inlet and outlets)
- Ecotoxicity of special waste and evaluation of the ecotoxicological risk of polluted sites
- Control of leachate from landfills, incineration residues, WWTP sludges, dredged sediments, ...

Microbiotests Toxkits


FRESHWATER / WASTEWATER

Tests with crustaceans



DAPHTOKKIT F

24-48h mobility inhibition test
with *Daphnia magna*
Conform with ISO Standard 6341 and
OECD Guideline 202



THAMNOTOKKIT F

24h mortality test
with *Thamnocephalus platyurus*
Conform with ISO Standard 14380



CERIODAPHTOKKIT F

24h mortality test
with *Ceriodaphnia dubia*
Adheres to test method US EPA 2002.0




Microbiotests Toxkits

Tests with micro-algae

ALGALTOXKIT F

72h growth inhibition test
with *Pseudokirchneriella subcapitata*
Conform with ISO Standard 8692 and
OECD Guideline 201



Tests with higher aquatic plants

DUCKWEED TOXKIT F

72h growth inhibition test
with *Spirodela polyrhiza*
Conform with ISO Standard 20227




Microbiotests Toxkits


Tests with rotifers



24h mortality test
with *Brachionus calyciflorus*
Conform with ISO Standard 19827



48h reproduction test
with *Brachionus calyciflorus*
Conform with ISO Standard 20666



Tests with ciliated protozoans



24h growth inhibition test
with *Tetrahymena thermophila*
Standardised screening method




Microbiotests Toxkits

SEAWATER / ESTUARIES

Tests with micro-algae



72h growth inhibition test
with *Phaeodactylum tricornutum*
Conform with ISO Standard 10253



Tests with crustaceans




24h mortality test
with *Artemia franciscana*
In current use in several countries as a
recommended acute toxicity test



Tests with rotifers



24h mortality test
with *Brachionus plicatilis*
Conform with ISO Standard 19820



Microbiotests Toxkits

SEDIMENTS / SOILS / SOLID WASTES

Tests with crustaceans



OSTRACODTOXKIT F

6 days mortality and growth inhibition test
with *Heterocypris incongruens*
Conform with ISO Standard 14371



Tests with terrestrial plants



PHYTOTOXKIT SOLID SAMPLES

PHYTOTOXKIT LIQUID SAMPLES



72h germination and root growth inhibition test
with seeds of 3 higher plants
Conform with ISO Standard 18763





ISO 11348-3: Water quality

Determination of the inhibitory effect of water samples on the light emission of *Aliivibrio fischeri* (luminescent bacteria test)




Freeze-dried reagent of *Aliivibrio fischeri*
(strain NRRL B-11177)



ISO 8692 : Water quality

Fresh water algal growth inhibition test with unicellular green algae




Microbiotests ALGALTOXKIT F

Standardized and culture-independent
toxicity testing and research :

Using algae obtained from alginate beads :
Pseudokirchneriella subcapitata /
Desmodesmus subspicatus






ISO 6341 : Water quality

Determination of the inhibition of mobility of *Daphnia magna* Straus (Cladocera, Crustacea) - Acute Toxicity Test



*Standardized and culture-independent
toxicity testing and research :*

*Using *Daphnia* hatched directly from ephippia*



The
New
Generation In
Acute Toxicity
Testing

AQUA
SCIENCE

Water Quality Luminescent Bacteria Test

determination of the inhibitory effect of water samples on the light emission of *Aliivibrio fischeri*

ISO 11348 consists of the following parts :


Part 1 : Method using freshly prepared bacteria (ISO 11348-1)

Part 2 : Method using liquid dried bacteria (ISO 11348-2)

Part 3 : Method using freeze dried bacteria (ISO 11348-3)



Aliivibrio fischeri (formerly *Vibrio fischeri*) is a marine bacterium that luminesces as a natural part of its metabolism under optimal environmental conditions.



Principle Luminescent bacteria test

When exposed to a toxic substance, the metabolic respiratory process of this bioluminescent bacterium is disrupted, reducing light output.

This inhibition of luminescence by *Aliivibrio fischeri* can be measured after 5-30 minutes with a luminometer.

The reduction of light intensity measured, directly correlates with the degree of toxicity of the sample relative to the control sample.



TOXIN CONCENTRATION

FIELDS OF APPLICATION

LUMINESCENT BACTERIA TESTS

☐ **Environmental analysis**

- Monitoring wastewater
- Screening sediment quality and soils
- Hazardous waste classification
- Bioremediation
- Surveillance environmental emergencies

☐ **Drinking water production**

- Detection of deliberate or accidental contamination
- Control of water supply sources (ground- and surface water)
- Checks of the distribution network and delivery systems

☐ **Waste Water Treatment Plants**

- Screening of influent before treatment
- Cost containment due to overload of toxic substances
- Protecting biological treatment
- Controlling the cost of chemical additives
- Effluent discharge legislation

☐ **Industrial applications**

- Drilling fluids / Petrochemical industry
- Mining wastes, soil and water
- Industrial effluents

Biosensor test using bioluminescent bacteria have been used for more than 25 years and their capability in detecting toxic substances is well understood.



DISCONTINUED



Highest Quality Bioluminescent Reagent



**Freeze-dried reagent with *Aliivibrio fischeri*
(conform ISO 11348-3)**




One hundred million cells / vial

Store frozen at -20 °C to – 25°C

Useful up to 4 hours after reconstitution

Two year shelf life from date of manufacture

Sensitive to > 3.600 chemical compounds



CONSUMABLES

BioLight Recon

BioLight Diluent

BioLight Solid Phase

Diluent

BioLight Salinity

Adjustment





VERSATILE combined instrument benchtop and field reader

EXPANDABLE can accommodate up to 2 cooling blocks


COMPLETE on-board instruction software and data evaluation

CUSTOMIZABLE instructions and protocols

SPEED reading time only 3 seconds

ACCESSIBLE wireless instrument control







BioLight Cooling blocks

Liquid (12 mm tubes) and Solid Tests (17 mm tubes)



Sample well
 $15 \pm 0,5 \text{ }^{\circ}\text{C}$

Possibility of working with two
cooling blocks.
(Up to 60 wells at the same time)

Custom PLC Software Interface on board

Easy to Set Up

Built-in test methods and protocols

Customizable instructions and settings

Multiple languages

Remote access and assistance

Free online software update up to two years after purchase

EU service center in co-operation with Ecotox Lds.



BioLight Test Results

Protocol:	Quality Test
Method:	Zinc
BioLight Lot #:	AS10970422
User:	admin
Description:	2022/07/26/16:16-BL10970422,Zinc,Run1,5P

Number of Controls:	1
Number of Dilutions:	4
Number of Samples:	1
Duplicates:	2
Incubation 1 (min):	15

	1	2	3	4	5
A	C	D	D	D	S
B	R	R	R	R	R
C	R	R	R	R	R
D					
E					
F					


AQUA SCIENCE [Back](#)



ACUTE TOXICITY TEST PROTOCOLS

- **EC50 Test - 2%, 45% or 81,9%**
Testing that provides a calculated EC50 value an measure relative toxicity
- **% Effect Test - 2%, 45% or 81,9%**
Testing using a single concentration
- **Equivalence Test**
Used with an undetermined sample to compare the relative acute toxicity with a control sample
- **ISO**
Complies with standardized method ISO11348-3
- **Solids Test**
Used for soil and sediments
- **Phenol Standard**
- **Zinc Standard**
- **3,5 – Dichlorophenol Standard**





AQUA
SCIENCE

Interested in more information?

www.microbiotests.com